


## STANDARD DIMENSIONS

PART NUMBER	Bottle Material	Closure Material	OVERFLOW (BRIM) CAPACITY (mL)	A IN (MM)	B IN (MM)	Ø C IN (MM)
160-04-0125	PETG	HDPE	175	4.28" (108.7)	2.06" (52.3)	1.72" (43.7)
160-04-0125-S	PETG	HDPE	175	4.28" (108.7)	2.06" (52.3)	1.72" (43.7)
160-04-0250	PETG	HDPE	320	5.79" (147.1)	2.31" (58.7)	1.72" (43.7)
160-04-0250-S	PETG	HDPE	320	5.79" (147.1)	2.31" (58.7)	1.72" (43.7)
160-04-0500	PETG	HDPE	610	7.05" (179.1)	2.88" (73.2)	1.72" (43.7)
160-04-0500-S	PETG	HDPE	610	7.05" (179.1)	2.88" (73.2)	1.72" (43.7)
160-04-1000	PETG	HDPE	1200	8.75" (222.3)	3.62" (91.9)	1.72" (43.7)
160-04-1000-S	PETG	HDPE	1200	8.75" (222.3)	3.62" (91.9)	1.72" (43.7)



REV.  
A

TITLE:  
PURILLEX® PETG BOTTLE  
DIMENSIONS