


**STANDARD DIMENSIONS**

<b>PART NUMBER</b>	<b>VOLUME</b>	<b>A</b> <b>INCHES (MM)</b>	<b>∅ B</b> <b>INCHES (MM)</b>	<b>C</b> <b>INCHES (MM)</b>	<b>D</b> <b>INCHES (MM)</b>
271-01-905-50	.5 ml	0.65" (16.5)	0.31" (7.8)	0.08" (2.0)	1.22" (30.9)
271-01-915-50	1.5 ml	0.71" (17.9)	0.42" (10.7)	0.09" (2.2)	1.65" (41.8)



TITLE:  
Purillex™ Microcentrifuge  
Vial Web Drawing

Rev.  
A