

STANDARD DIMENSIONS						
	PART NUMBER	NOMINAL VOLUME (mL)	OVERFLOW (BRIM) CAPACITY (mL)	A IN (MM)	B IN (MM)	ϕ C IN (MM)
FIGURE 1	150-01-0100	100	165	4.27" (108.5)	2.12" (53.8)	2.19" (55.6)
	150-01-0250	250	380	5.75" (146.1)	2.12" (53.8)	2.75" (69.9)
	150-01-0500	500	545	6.33" (160.9)	2.12" (53.8)	3.35" (85.1)
	150-01-1000	1000	1200	8.60" (218.4)	2.12" (53.8)	3.98" (101.1)
	150-01-2000	2000	2168	9.77" (248.1)	2.12" (53.8)	5.01" (127.3)
FIGURE 2	150-01-0050	50	62	3.31" (84.2)	1.52" (38.7)	1.44" (36.6)

FIGURE 1





FIGURE 2



TITLE: **LAB BOTTLE DIMENSIONS**

Rev. D