



**STANDARD DIMENSIONS**

	<b>PART NUMBER</b>	<b>Material</b>	<b>NOMINAL VOLUME (mL)</b>	<b>OVERFLOW (BRIM) CAPACITY (mL)</b>	<b>A IN (MM)</b>	<b>B IN (MM)</b>	<b>Ø C IN (MM)</b>
<b>FIGURE 1</b>	170-01-0100	PFA	100	165	4.27" (108.5)	2.12" (53.8)	2.19" (55.6)
	170-01-0100-A	PFA	100	165	4.27" (108.5)	2.12" (53.8)	2.19" (55.6)
	170-01-0250	PFA	250	380	5.75" (146.1)	2.12" (53.8)	2.75" (69.9)
	170-01-0250-A	PFA	250	380	5.75" (146.1)	2.12" (53.8)	2.75" (69.9)
	170-01-0500	PFA	500	545	6.33" (160.9)	2.12" (53.8)	3.35" (85.1)
	170-01-0500-A	PFA	500	545	6.33" (160.9)	2.12" (53.8)	3.35" (85.1)
	170-01-1000	PFA	1000	1200	8.60" (218.4)	2.12" (53.8)	3.98" (101.1)
	170-01-1000-A	PFA	1000	1200	8.60" (218.4)	2.12" (53.8)	3.98" (101.1)
	170-01-2000	PFA	2000	2240	9.77" (248.1)	2.12" (53.8)	5.01" (127.3)
	170-01-2000-A	PFA	2000	2240	9.77" (248.1)	2.12" (53.8)	5.01" (127.3)
	170-02-0100	FEP	100	162	4.27" (108.5)	2.12" (53.8)	2.19" (55.6)
	170-02-0250	FEP	250	352	5.75" (146.1)	2.12" (53.8)	2.75" (69.9)
	170-02-0500	FEP	500	588	6.33" (160.9)	2.12" (53.8)	3.35" (85.1)
	170-02-1000	FEP	1000	1182	8.60" (218.4)	2.12" (53.8)	3.98" (101.1)
	170-02-2000	FEP	2000	2240	9.73" (247.0)	2.12" (53.8)	5.05" (127.3)
	170-03-0250-G	ETFE	250	380	5.86" (148.8)	2.12" (53.8)	2.75" (69.9)
170-03-0500-G	ETFE	500	545	6.44" (163.5)	2.12" (53.8)	3.35" (85.1)	
170-03-1000-G	ETFE	1000	1200	8.71" (221.2)	2.12" (53.8)	3.98" (101.1)	
<b>FIGURE 2</b>	170-01-0050	PFA	50	58	3.31" (84.2)	1.52" (38.7)	1.44" (36.6)
	170-01-0050-A	PFA	50	58	3.31" (84.2)	1.52" (38.7)	1.44" (36.6)
	170-02-0050	FEP	50	62	3.31" (84.2)	1.52" (38.7)	1.44" (36.6)



**FIGURE 1**



**FIGURE 2**



TITLE: **PURILLEX BOTTLE DIMENSIONS**