



### STANDARD DIMENSIONS

PART NUMBER	A INCHES (MM)	$\phi$ B INCHES (MM)	$\phi$ C INCHES (MM)
550-740-400891	2.79" (70.8)	1.00" (25.4)	6.00" (152.4)
550-740-400892	1.64" (41.75)	1.00" (25.4)	6.00" (152.4)
550-740-400893	1.53" (38.8)	1.00" (25.4)	6.00" (152.4)
550-740-400894	0.84" (21.3)	1.00" (25.4)	6.00" (152.4)
550-740-400895	0.75" (19.1)	1.00" (25.4)	6.00" (152.4)
550-740-400896	0.63" (15.9)	1.00" (25.4)	6.00" (152.4)




 10321 West 70th St.  
 Eden Prairie, MN  
 55344-3446 USA  
 (952)935-4100

TITLE:  
**VIAL CLEANING SYSTEM,  
 RACK DIMENSIONS**